

## SEM : Quanta FEI

- Log into computer
- Open xT microscope server
  - Wait for the three buttons to turn green
  - Hit “Start”
- Log-into microscope with same log-in as the computer
  - Home stage
  - Note any error message and clear message screen
  - Minimize server window by double clicking in blue bar




### General Notes:

- Vacuum Status:



Status	
Chamber Pressure:	1.45e-4 Pa
Gun Pressure:	2.67e-7 Pa
Emission Current:	152 µA

**Green:** PUMPED to the desired vacuum mode  
**Orange:** TRANSITION between two vacuum modes  
(pumping / venting / purging)  
**Grey:** VENTED for sample or detector exchange

-   
- Short cut keys:
  - +/- = change magnification
  - \* = rounds off magnification to even number
  - F2 = takes picture
  - F3 = videoscope
  - F5 = maximizes/minimizes window
  - F7 = reduced area focus
  - Left click = focus
  - Left click + shift = stigmat
  - ESC = cancel button, will stop whatever the active command is
- Saving and retrieving files
  - Save all files on the left computer in your user account
  - Do not put anything in the right, microscope computer as this will void the warranty
- On first start up – go to tools → preferences → general → lower z upon vent  
→ data bar → date → move to top


### General order of operations for the SEM

1. Insert sample/pump down
2. Rough focus samples
  - Link to z
  - Align microscope
3. Map samples
4. Fine focus and image
  - Fine stigmata each image
5. Return all settings to normal and remove sample

### SEM Operation


## 1. Load sample





- always wear gloves when touching anything to go into the chamber
- “Vent”
  - “Yes” – you want to vent
- Insert sample, and screw down
- Hold door shut and “Pump”
- Wait for icon to turn completely green
- Change the sample height
  - Press and hold the scroll button and move the mouse up to move the stage to ~10 mm. (for a flat sample on the standard sample 25 mm stub = 55 mm on z)
  - After your initial focus and z-height alignment, it can be moved closer



## 2. Rough Focus and Alignment

- Good conditions to start = 10 kV and spot size = 1-2
- Click on window that you want to view samples (bar will turn blue) and unpause (||)
- “Beam On”








- Rough focus the image (should be around WD = 10mm)
- Link Z (  ) → 
  - Tells the instrument where the stage height actually is.
  - Re-link Z often during imaging 
- Lens Align 
  - This should be done at a mag greater than your normal working mag
  - Want this to blink in and out, not move up and down/side to side
- Crossover
  - Spot size should be centered over the X
  - We do not have control over this feature, but it should be checked after changing spot size or kV



## 3. Map Samples

## 4. Fine focus and Imaging

- After focusing sample and re-linking z, can move sample closer to gun for higher res imaging
- Sigmation = shift and right click
  - “Focus” removes the direction
  - “Stigmat” sharpens image

- Think of focus and stigmata as working together: focus then stigmatate then refocus and restigmatate until you are happy with the quality of the image
- Generally, you shouldn't have to make significant changes in the stigmatation
- To use reduced area screen use: 
- To move image, press and hold scroll button and move mouse the direction you want to go
- Double click will center object that you double clicked
- Scroll Brightness and Contrast buttons to adjust/use videoscope () to aid
- Use reduced area to fine focus and stigmata 
- To take picture, hit  or F2 and save files when prompted in your folder
- Use () to change the scan speed or the arrows

#### 5. Return all settings to normal and remove sample

- “Beam On” = turn beam off (wait for valve release)
- Vent
- Remove sample and any accessory, return to high vacuum mode if necessary
- “Pump” – wait for () to turn green
- Log out of microscope screen
- File → Exit out of microscope screen
- If no one else is coming:
  - i. Turn off server → “Stop” → wait for blue bar to go away
  - ii. Exit out of server (right click on blue bar)
  - iii. Log off computer
- **Low vac mode:**
  - Vent
  - Install low vac accessory on the end of gun
    - it will gently push in like a button
    - use CCD to verify that it is pushed in evenly
  - Change to “low vac”
  - Change chamber pressure to 0.6 to 0.68 Torr to start
  - Pump
  - Choose “low vac” accessory option
  - Vacuum will pump, then purge, then pump again before () turns green
  - Can change the vacuum pressure during imaging to optimize
  - Start at 8 kV

## Sputter Coater


1. Turn on gas tank
2. Turn on sputter coater (switch on back)
3. → System start
4. → Screens
5. → Manual Sputter set up
  - a. → Turn on mechanical pump (flip valve on pump to make sure it is on)  
Wait until pressure is down to 40 mTorr
  - b. → Open gas valve  
Adjust gas pressure to 50-70 mTorr using valve
6. → Set point  
Should be 35-40
7. Turn on thickness monitor
8. → Rotation on
9. → Power on  
Let run until desired thickness reached (no more than 10 minutes)
10. → Power off
11. → close gas valve
12. → turn off mechanical pump  
Turn valve to closed
13. Turn off gas
14. Turn off sputter coater
15. Turn off thickness monitor

## **Oxford INCA EDX**

\*\*\*Use uncoated Samples for EDX

1. Start the SEM as normal
  - 10 mm is the optimum working distance for the geometry of the detector during EDX measurements
  - Use a KV of 2x that of the highest energy you hope to see (found on slide rule or mouse pad by the computer)
  - Start with a spot size of ~4
2. Focus and stigmata on your sample at this conditions
3. Pause the CCD camera
4. Open the INCA program on the support computer
  - INCA can be operated under three different areas. Each area has its own flow chart to follow during analysis
    - (1) Analyzer – gathers information from the whole field of view
    - (2) Point and ID – can choose where you want to get information
    - (3) Mapping – can map the spectra
5. Image the area of interest
6. Log off SEM as usual and close INCA program


## \*\* Analyzer \*\*



- **Project**
  - Name project or start a new project by choosing the thumb tack icon
- **Sample**
  - Name your sample or start a new sample by choosing the thumb tack icon
  - Each project can have multiple samples
- **Microscope Set-Up**
  - Run the cycler
  - Input rate > 2k. This can be adjusted by increasing the spot size
  - Dead time should be <40%
  - Choose a process time between 1-6. 6 is the longest process time and will give the highest resolution; 1 is the shortest and will have the lowest resolution. 4-5 is a good place to start. This can be adjusted to influence the dead time
- **Acquire Spectra**
  - Each sample can have multiple points of interest
  - Choose an appropriate acquisition time and acquisition time and kV range

- **Confirm Elements**
  - Confirm the elements that may be present in your Sample
- **Quant**
  - Identify the correct elements
    - Quant set-up
      - "Transfer" the elements from your spectra or choose the element of interest from the drop down menu
- **Report**
  - Choose report format you want, export, and save the report as a word document
  - Repeat for each spectra

## \*\*Point and ID \*\*



- Project
  - Name project or start a new project by choosing the thumb tack icon
- Sample
  - Name your sample or start a new sample by choosing the thumb tack icon
  - Each project can have multiple samples
- Microscope Set-Up
  - Run the cyclotron
  - Input rate > 2k. This can be adjusted by increasing the spot size
  - Dead time should be <40%
  - Choose a process time between 1-6. 6 is the longest process time and will give the highest resolution; 1 is the shortest and will have the lowest resolution. 4-5 is a good place to start. This can be adjusted to influence the dead time
- Site of Interest
  - Image Set-up
    - Adjust resolution
  - Focus on your sample area under the microscope
  - Record image on INCA (green circle)

- Acquire Spectra
  - Each sample can have multiple points of interest
  - Choose an appropriate acquisition time and acquisition time and kV range
- Confirm Elements
  - Confirm the elements that may be present in your sample
- Quant
  - Identify the correct elements
  - Quant set-up

- “Transfer” the elements from your spectra or choose the element of interest from the drop down menu

➔ Report

- Choose report format you want, export, and save the report as a word document
- Repeat for each spectra